

# ANATOMY

## Artificial Bone Models

Nature is our Model



### QS 10/14 • ARTIFICIAL HUMAN SKELETON

Male, modelled according to nature, in SOMSO-PLAST<sup>®</sup>. Lifelike representation of all anatomical details of the bone structure. Skull with removable calvarium and mandible. Joints are flexibly mounted, the upper and lower extremities can be removed. The right and the left foot can be detached from the lower legs. With flexible vertebral column and spinal cord with nerve endings. Mounted standing upright with a stand for hanging with castors. Height: 180 cm (skeleton 171 cm), width: 55 cm, depth: 55 cm, weight: 11.5 kg