



## QS 10/8 • ARTIFICIAL HUMAN SKELETON

Female, modelled according to nature, in SOMSO-PLAST®. Lifelike representation of all anatomical details of the bone structure. Skull with removable calvarium and mandible. Joints are flexibly mounted, the upper and lower extremities can be removed. The right and the left foot can be detached from the lower legs. Mounted upright on a stand with castors. With a dust cover.

Height: 181 cm (skeleton 170 cm), width: 55 cm, depth: 55 cm, weight: 11 kg